



Issued to:

Anhui Intco Medical products Co.,Ltd  
No. 6, Haitang South Road  
Suixi Wuhu Modern Industrial Park  
Suixi County, Anhui Province  
Huaibei City  
235000  
China

Notified Body: 2777

SATRA customer number: P21130

## EU Type-Examination Certificate

**Certificate number: 2777/14815-04/E00-00**

This EU Type-Examination Certificate covers the following product group(s) supported by testing to the relevant standards/technical specifications and examination of the technical file documentation. It has been issued Under Module B of Regulation 2016/425 on personal protective equipment. This product group has been shown to satisfy the applicable essential health and safety requirements as a Category III product.

### Product reference:

Blue :  
NGPF7000(XS)-7005(XXL)  
NGV/B/H/P/XE100&200 13-18  
NGV/B/H/P/XE10L&20L 13-18  
SNV/B/H/PE100&200 13-18  
SNV/B/H/PE10L&20L 13-18  
Black:  
BNPF7000(XS)- 7005(XXL)  
NGV/B/H/P/XE100&200 43-48  
NGV/B/H/P/XE10L&20L 43-48  
SNV/B/H/PE100&200 43-48  
SNV/B/H/PE10L&20L 43-48  
Red:  
NGV/B/H/P/XE100&200 73-78  
NGV/B/H/P/XE10L&20L 73-78  
SNV/B/H/PE100&200 73-78  
SNV/B/H/PE10L&20L 73-78

### Description:

Disposable Nitrile Gloves/ SYNGUARD @Nirile Exam Gloves, Powder free  
Weight:  $\geq 3.0g$ , Thickness:  $\geq 0.06mm$

White :  
WNPF7000(XS)-7005(XXL)  
NGV/B/H/P/XE100&200 23-28  
NGV/B/H/P/XE10L&20L 23-28  
SNV/B/H/PE100&200 23-28  
SNV/B/H/PE10L&20L 23-28  
Pink:  
NGV/B/H/P/XE100&200 53-58  
NGV/B/H/P/XE10L&20L 53-58  
SNV/B/H/PE100&200 53-58  
SNV/B/H/PE10L&20L 53-58  
Green:  
NGV/B/H/P/XE100&200 83-88  
NGV/B/H/P/XE10L&20L 83-88  
SNV/B/H/PE100&200 83-88  
SNV/B/H/PE10L&20L 83-88

Violet:  
NGPF7000-V1(XS)-7005-V1(XXL)  
NGV/B/H/P/XE100&200 33-38  
NGV/B/H/P/XE10L&20L 33-38  
SNV/B/H/PE100&200 33-38  
SNV/B/H/PE10L&20L 33-38  
Yellow:  
NGV/B/H/P/XE100&200 63-68  
NGV/B/H/P/XE10L&20L 63-68  
SNV/B/H/PE100&200 63-68  
SNV/B/H/PE10L&20L 63-68  
Orange:  
NGV/B/H/P/XE100&200 93-98  
NGV/B/H/P/XE10L&20L 93-98  
SNV/B/H/PE100&200 93-98  
SNV/B/H/PE10L&20L 93-98

### Sizes:

6-11(XS-XXL)

### Classification:

**EN ISO 374-1:2016+A1:2018/Type B**  
40% Sodium hydroxide (K)  
30% Hydrogen peroxide (P)  
37% Formaldehyde (T)

Level	EN ISO 374-4:2019 Degradation %
6	-68.1
2	30.5
5	9.5

### EN ISO 374-5:2016

Protection against Bacteria and Fungi  
Protection against Viruses

### Level

Pass  
Pass

### Standards/Technical specifications applied:

EN ISO 21420:2020; EN ISO 374-1:2016+A1:2018; EN ISO 374-5:2016; EN 421:2010 (Radioactive contamination only)

### Technical reports/Approval documents:

SATRA: CHT0296241/2012, CHM0298100/2020/EN/A, CHM0298100/2020/EN/B, CHT0338202/2239

SGS: CH:TX:1142011147, CH:TX:1142011145-1, CH:TX:1142011148

TUV: 7191234075-CHM20-02-TSL, 7191235025-EEC20-WBH\_CR1, 721652920

Signed on behalf of SATRA:

Kayleigh Aylward

Date first issued: 20/07/2020

Date of issue: 25/04/2024

Expiry date: 20/07/2025

# TERMS AND CONDITIONS

The following conditions apply in addition to SATRA's standard terms and conditions of business and those given in the current certification agreement. This certificate has been issued in accordance with Annex V (Module B) of the applicable legislation (see note 11).

Please note:

1. Where the product is classified as category III then CE or UKCA Marking of production is reliant on current compliance with module C2 or Module D of the applicable legislation (See note 11). (Except that specifically produced to fit an individual user).
2. Full details of the scope of the certification and product(s) certified are contained within the manufacturer's technical documentation.
3. Where a translation of this certificate exists, the English language version shall be considered as the authoritative text.
4. Certification is limited to production undertaken at the sites listed in the manufacturer's technical documentation.
5. Ongoing manufactured product shall be consistent with the product(s) certified and listed on this certificate and an EU or UKCA declaration of product conformity shall be made available in accordance with the applicable legislation (See note 11).
6. The Manufacturer shall inform SATRA of any changes to the certified product or technical documentation.
7. Where results obtained during type testing are within the budget of uncertainty when compared to the pass requirement, classification, or performance level, then it is the responsibility of the manufacturer to ensure that the factory production control and manufacturing tolerances are such that the product placed on the market meets with the stated requirements, classifications or performance levels.
8. This certificate shall be kept together with the relevant technical documentation in a safe place by the client named on this certificate. Production of this certificate and other documentation may be required by a representative of the EC member state, or UK government.
9. This certificate relates only to the condition of the testable items at the time of the certification procedure and is subject to the expiry date shown.
10. SATRA reserves the right to withdraw this certificate if it is found that a condition of manufacture, design, materials, or packaging have been changed and therefore no longer comply with the requirements of the applicable legislation (See note 11).
11. These terms and conditions shall apply to the requirements set out in Regulation (EU) 2016/425 of the European Parliament and of the council of 9th March 2016 on personal protective equipment or to UK legislation relating to UKCA Marking as defined within the issued certificate.